

### Hardness/Iron Test Kit - \$55

Iron and Hardness

### Kit includes:

1 x 250 mL Plastic Bottle

### **Sampling Instructions:**

#### READ ALL INSTRUCTIONS BEFORE PROCEEDING WITH SAMPLE COLLECTION

\*Sampling video available on our website at www.nlslab.com under support material

Collect water that precedes any water treatment within your home. This is representative of the well water itself.

- 1. Collect water from selected location, avoiding new plumbing if possible.
- 2. Remove aeration screen from faucet if one is present.
- 3. Run water at maximum flow for at least 2 minutes.
- 4. Fill the containers.

# Follow the return shipping instructions to ensure kit is shipped back to the lab properly:

- Fill out the required information below, and label your sample bottles with date/time of collection.
- Samples are accepted 7:00 am to 4:00 pm Monday Friday at the Crandon lab and 8:00 am to 4:00 pm Monday Friday at the Waukesha lab. NLS receives deliveries from UPS, Fed Ex, SpeeDee Delivery, USPS, and other private courier services. Use any of these courier services to return samples to NLS.
- If you have any questions regarding this sampling kit, please call our Client Service Department at 800-278-1254

## **Hardness Panel Test Kit - \$55**

Iron and Hardness

## Return your sample as soon as possible to one of our two locations:

Northern Lake Service
400 N Lake Ave
Crandon, WI 54520

Northern Lake Service
2420 N Grandview Blvd
Waukesha, WI 53188

## Please provide the following information:

Name:	
Address:	
City/State/Zip:	
Phone:	Email:
Sample Collection Date:	Sample Collection Time:AM/PM
Sample Collection Location (ex. Pressure Ta	nk)
Sample Collected By:	
Report by email only Report by n	nail only Report by email and mail