

Nitrate Test Kit - \$31.50

Kit includes:

- 1 x 1 Ziplock bag TO BE FILLED WITH ICE PRIOR TO SHIPMENT
- 1 x 125 mL Plastic Bottle

General Sampling Instructions: (Choose 1 of 2 options)

READ ALL INSTRUCTIONS BEFORE PROCEEDING WITH SAMPLE COLLECTION *Sampling video available on our website at www.nlslab.com under support material

Option 1: Collect water that is representative of the way you actually consume your water.

Option 2: Collect water that precedes any water treatment within your home. This is representative of the well water itself.

- 1. Collect water from selected location, avoiding new plumbing if possible.
- 2. Remove aeration screen from faucet if one is present.
- 3. Run water at maximum flow for at least 2 minutes
- 4. Fill the container.

Follow the return shipping instructions to ensure kit is shipped back to the lab properly:

- Fill out the required information below, and label your sample bottles with date/time of collection.
- Samples are accepted 7:00 am to 4:00 pm Monday Friday at the Crandon lab and 8:00 am to 4:00 pm Monday Friday at the Waukesha lab. If you are shipping your kit to NLS, samples MUST arrive at the laboratory no longer than 28 hours from the time of collection. Do not ship on Friday.
- NLS receives deliveries from UPS, Fed Ex, SpeeDee Delivery, USPS, and other private courier services. Use any of these courier services to return sample to NLS.
- If you have any questions regarding this sampling kit, please call our Client Service Department at 800-278-1254

Nitrate Test Kit - \$31.50

Return your sample as soon as possible to one of our two locations:

Northern Lake Service
400 N Lake Ave
OR
Crandon, WI 54520

Northern Lake Service
2420 N Grandview Blvd
Waukesha, WI 53188

Please provide the following information:

Name:		
Address:		
City/State/Zip:		
Phone:		
Sample Collection Date:	Sample Collection Time:	AM/PM
Sample Collection Location (ex. Pressu	re Tank)	
Sample Collected By:		
Report by email only Repor	t by mail only Report by em	nail and mail